

Software collaboration

ETH + QuTech

19-5-2016

Thursday

Goal: Understand requirements, capabilities and potential for collaboration.

8:30 - 9:00 arrival at QuTech + Coffee

9:00 - 10:00 Presentations by Johannes and Adriaan

Current/desired framework Johannes Heinsoo

Current Delft framework (PycQED + QCodes)

11:00 - 12:00 Discussion to

1. Summarize conclusions of initial presentations

2. Identify sub-goals for rest of the visit

12:00- 13:00 lunch

13:00 - 15:00 Talk with other lab-members on data acquisition software (Ramiro and Chris in particular)

17:00 - 18:00 Meeting Johannes and Adriaan, draw conclusions from the day.

19:00 - Dinner, opportunity to discuss software in a more relaxed/less structured setting

Friday

Goal: plan for software collaboration in QuSurf + decision on platform

9:00 - 10:00 Meeting Johannes + Adriaan, goal: a plan for structural software collaboration

12:00 - 12:30 Werkbepreking lunch

12:30 - 13:30 Werkbepreking

16:00 Meeting with Leo, Johannes, Niels and Adriaan to discuss next steps on software platform. Additionally present plan for collaboration.

17:00 Celebrate @ TPKV!

Goal of this presentation

Understand

- **Requirements** of data acquisition software
- **Current capabilities** and platform
- Potential for **collaboration**


Requirements

- Controls the setup, logs and plots the data
- Automatable and extensible
- Easy to learn/use
- Grows with newest experimental insights



Current Platform: Building layers of abstraction

- High level - 'Quantum Compiler'
 - Solid?
- Mid level - System specific bridge to abstraction
 - PycQED
- Low level - Instrument control software
 - QCodes/QTlab



The current framework is build around a minimal set of concepts

Core concepts

- Parameter
 - Sweep function
 - Detector function
- Instrument
- MC/Loop
- Analysis

Derived concepts

- Composite parameter
- hard/soft sweeps/parameters
- Meta-instrument
- Adaptive measurements
- Qubit-object

Every experiment consists of a Loop

Example: Heterodyne experiment

Loop:

1. Some **parameter**(s) is/are varied
2. Some **parameter** is measured
3. Data is saved and analyzed

PycQED syntax

```
1 MC.set_sweep_function(source.frequency)
2 MC.set_sweep_points(np.linspace(start, stop, steps))
3 MC.set_detector_function(HeterodyneDetector())
4 MC.run(name='Heterodyne')
5 ma.MeasurementAnalysis()
```

QCodes syntax

```
7 data = qc.Loop(source.frequency[start:stop:step]).each(
8     ... Heterodyne.IQ).run(name='Heterodyne')
9 qc.QTPlot(data.IQ)
```

The Loop/MC takes care of
standardised datasaving, logging
and live plotting

The notion of a parameter is sufficiently abstract that it allows nested/composite measurements



Example: T1 measurement

```
MC.set_sweep_function(swf.Flux_Control_mV(dac_channel))
MC.set_sweep_points(virtual_dac_voltages)
MC.set_detector_function(cdet.T1_Detector(qubit))
MC.run(measurement_name='T1_Mux0_Bottom_Qubit',
      debug_mode=False)]
```

Consists of :

1. Finding the resonator
2. Finding the qubit
3. Calibrating the pulse amplitude
4. Performing a T1 measurement

A **composite parameter** can contain other parameters and/or complete loops with analysis




The concept of the loop is sufficiently flexible that it allows adaptive measurements

Example: Numerically optimised mixer skewness

```
ad_func_pars = {'adaptive_function': nelder_mead,  
                'x0': [1.0, 0.0],  
                'initial_step': [.15, 10],  
                'no_improv_break': 5,  
                'minimize': True,  
                'maxiter': 500}  
MC.set_sweep_functions([S1, S2])  
MC.set_detector_function(d)  
MC.set_adaptive_function_parameters(ad_func_pars)  
MC.run(name=name, mode='adaptive')  
a = MA.OptimizationAnalysis()
```

The adaptive loop can use **any adaptive function** on the results of anything we can quantify/measure




Meta-instruments allow layers of abstraction


Example: Flux-control
Converts fluxes to dac-voltages
using a calibrated correction matrix

Example: Composite AWG
Acts as a single multi-channel AWG
but talks to underlying instruments


A **meta-instrument** can
contain other instruments
but acts like a regular
instrument



The **qubit object** is a special meta-instrument that executes small scripts (e.g. find frequency using Ramsey) and holds parameters



Potential for collaboration



Software collaboration

ETH + QuTech

19-5-2016

AWG sequencing (made easy)

Pulses(dict) containing pulses by
“name” (X180 etc)

Pulse_pars (dict)

AWG segment is made by adding
pulses (by key) in a for loop

Pulse-definition is contained in
pulse-lib (currently only RO
pulse and SSB Drag pulse)

- Includes fix-point correction
- Automated sideband modulation
- Allows for pulse definitions including markers
- In-principle extensible to arbitrary nr of channels

```
MC.set_sweep_function(awg_swf.Randomized_Benchmarking(
    pulse_pars=self.pulse_pars, RO_pars=self.RO_pars,
    nr_cliffords, nr_seeds=nr_seeds))

pulses = {'I': deep,
          'X180': de,
          'mX180': d,
          'X': 149,
          'm': 150,
          'Y': 151,
          'm': 152,
          'Y': 153,
          'm': 154,
          'Y': 155,
          'm': 156,
          'Y': 157,
          'm': 158,
          'Y': 159,
          'm': 160}

def get_pulse_pars(self):
    self.pulse_pars = {
        'I_channel': self.pulse_I_channel.get(),
        'Q_channel': self.pulse_Q_channel.get(),
        'amplitude': self.amp180.get(),
        'sigma': self.gauss_sigma.get(),
        'nr_sigma': 4,
        'motzoi': self.motzoi.get(),
        'mod_frequency': self.f_pulse_mod.get(),
        'pulse_separation': self.pulse_separation.get(),
        'phase': 0,
        'pulse_type': 'SSB_DRAG_pulse'}
```

```
else:
    cl_seq = rb.randomized_benchmarking_sequence(n_cl)
    pulse_keys = rb.decompose_clifford_seq(cl_seq)
    pulse_list = [pulses[x] for x in pulse_keys]
    pulse_list += [RO_pars]
    el = multi_pulse_elt(i, station, pulse_list)
    el_list.append(el)
    seq.append_element(el, trigger_wait=True)
```